

## Product Datasheet

# Multi-species IL-5 (Bovine, Caprine, Ovine) Recombinant Protein (orb1820834)

<b>Description</b>	The Bovine/Caprine/Ovine IL-5 yeast-derived recombinant protein is not tagged and is naturally endot...
<b>Target</b>	IL-5
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VESTMNRLVA ETLTLLSSHR TLLIGDGNLM IPTPQHTNHQ LCIEEVFQGI DTLKNQTAQG DAVKKIFQNL SLIKEYIDLQ KRKCGGERWR VKQFLDYLQV FLGVINTEWT MES (113)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	13.0 kDa
<b>Entrez</b>	<a href="#">280825</a>
<b>Expiration Date</b>	6 months from date of receipt.