

Product Datasheet

Swine IL-6 Recombinant Protein (orb1820835)

Description	The Swine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GRLEEDAKGDATSDKMLFTSPDKTEELIKYILGKISAMRKEMCEKYKENSKEVL AENNLNLPKMAEKDGCFSQGFNQETCLMRITTTGLVEFQIYLDYLQKEYESNKG VEAVQISTKALIQTLRQKGNPKATTPNPTTNAAGLLDKLQSQNEWMKNTKIILIL RSLEDFLQFSLRAIRIM (183)
Purity	98%
Source	Yeast
MW	20.9 kDa
Entrez	399500
Expiration Date	6 months from date of receipt.