

Product Datasheet

Rabbit IFN gamma Recombinant Protein (orb1820840)

Description	The Rabbit IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QDTLTRETEH LKAYLKANTS DVANGGPLFL NILRNWKEES DNKIIQSQIV SFYFKLFDNL KDHEVIKKSM ESIKEDIFVK FFNSNLTKMD DFQNLTRISV DDRLVQRKAV SELSNVLNFL SPKSNLKKRK RSQTLFRGRR ASKY (144)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	100008602
Expiration Date	6 months from date of receipt.