

Product Datasheet

Swine IL-17A Recombinant Protein (orb1820847)

Description	The Swine IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIMIPQSPGCPKTEDKNFPQHVRVNLNILNRSTPARRPSDYSKRSTSPWTLQRNE DPERYPSVIWEAKCSHSGCINAEGKEDHHMNSVPIQQEILVLRREPRHCPNSFRL EKVMVTVGCTCVTPIVRHIS (130)
Purity	98%
Source	Yeast
MW	14.9 kDa
Entrez	449530
Expiration Date	6 months from date of receipt.