

Product Datasheet

Chicken TNFSF15 Recombinant Protein (orb1820855)

Description	The Chicken TNFSF15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	TNFSF15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ERSSHFLKQRAVA AVTDTLPSAEKPRAHLTVKKQEPSSTTGSHLPILQWEDKRGL AFTKNNLSYSSNALVIPVSGDYVYAQVTRFGPSDTSSKTSSVTAVITKVTDSTYPE PTQLLTSTKTLSEERNNWFQPIYLGAVVSLEIGDKLMVNVSDIKLVDYTKCHKTF GAFLL
Purity	98%
Source	Yeast
MW	19.1 kDa
Entrez	417247
Expiration Date	6 months from date of receipt.