

Product Datasheet

Chicken IL-4 Recombinant Protein (orb1820860)

Description	The Chicken IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LQLSVPLMES IRIVNDIQGE VSCVKMNVTD IFADNKTNNK TELLCKASTI VWESQHCHK LQGLFLNMRQ LLNASSTSLK APCPTAAGNT TSMEKFLADL RTFFHQLAKN K
Purity	98%
Source	Yeast
MW	12.4 kDa
Entrez	416330
Expiration Date	6 months from date of receipt.