

Product Datasheet

Swine IL-8 Recombinant Protein (orb1820862)

Description	The Swine IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARVSAELRCQ CINTHSTPFH PKFIKELRVI ESGPHCENSE IIVKLVNGKE VCLDPKEKWV QKVVQIFLKR TEKQQQQQ (78)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	396880
Expiration Date	6 months from date of receipt.