

## Product Datasheet

### Swine IL-8 Recombinant Protein (orb1820864)

<b>Description</b>	The Swine IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
<b>Target</b>	IL-8
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	ARVSAELRCQ CINTHSTPFH PKFIKELRVI ESGPHCENSE IIVKLVNGKE VCLDPKEKWV QKVVQIFLKR TEKQQQQQ (78)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	9.1 kDa
<b>Entrez</b>	<a href="#">396880</a>
<b>Expiration Date</b>	6 months from date of receipt.