

Product Datasheet

Equine IL-5 Recombinant Protein (orb1820869)

Description	The Equine IL-5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVESPMNRLVAETLTLSTHRTLLIGDGNLMIPTPEHKNHQLCIEEVFQIDTLKN QTVQGDVAKLFQNLSLIKGYIDLQKKKCGGERWRVKQFLDYLQEFLGVINTEW TIEG
Purity	98%
Source	Yeast
MW	12.9 kDa
Entrez	100034199
Expiration Date	6 months from date of receipt.