

Product Datasheet

Bovine IL-1 alpha Recombinant Protein (orb1820883)

Description	The Bovine IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SNVKYNFMRVIHQECILNDALNQSIIRDMSGPYLTATTLNMLEEAVKFDMVAYVSEEDSQLPVTLRISKTLQFVSAQNEDEPVLLKEMPETPKIIKDETLLFFWEKHGSM DYFKSVAHPKLFIA TKQEKLVMASGPPSITDFQILEK
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	281250
Expiration Date	6 months from date of receipt.