

Product Datasheet

Canine IL-1 beta Recombinant Protein (orb1820885)

Description	The Canine IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAMQSV DCKLQDISHKYLVL S NSYELRALHLNGENVNKQVVFHMSFVHGDESN NKIPVVLGIKQKNLYLSCVMKDGKPTLQLEKVDPKVYPKRKMEKRFVFNKIEIKNT VEFESSQYPNWIYISTSQVEGMPVFLGNTRGGQDITDFTMEFSS
Purity	98%
Source	Yeast
MW	17.5 kDa
Entrez	403974
Expiration Date	6 months from date of receipt.