

Product Datasheet

Swine TNF alpha Recombinant Protein (orb1820888)

Description	The Swine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SSSQTS DKPVAHVVANVKAEGQLQWQSGYANALLANGVKLKDNQLVVPTDGL YLIYSQVLFRRGQGPSTNVFLTHTISRIAVSYQTKVNLLSAIKSPCQRETPEGAEAK PWYEPIYLGGVFQLEKDDRLSAEINLPDYLDFAESGQVYFGIIAL
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	397086
Expiration Date	6 months from date of receipt.