

Product Datasheet

Equine IL-17A Recombinant Protein (orb1820893)

Description	The Equine IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIVIPQNPECPNTGDKNFPQNVKINLNLNRKTNSRRASDYHNRSTSPWNLHRN EDPERYPSVIWEAKCRHLGCVNAEGKVD FHMNSVPIQQEILVLRRESQNCPHSF QLEKMLVAVGCTCVTPIVRHMG
Purity	98%
Source	Yeast
MW	14.9 kDa
Entrez	100034142
Expiration Date	6 months from date of receipt.