

Product Datasheet

Bovine IL-12 p40 Recombinant Protein (orb1820894)

Description	The Bovine IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	IL-12p40
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MWELEKNVYVVELDWYPDAPGETVVLTCDTPEEDGITWTSQSDQSSEVLGSGKLT TIQVKEFGDAGQYTCHKGGEALSRSLLLLHKKEDGIWSTDILKDQKEPKAKSFLK CEAKDYSGHFTCWWLTAISTDLKFSVKSSRGSSDPRGVTCGAALLSAEKVSLEH REYNKYTVECQEGSACPAAEESLLIEVVVEAVHKLKYENYTSFFIRDIIKPDPPKN LQLRPLKNSRQVEVSWEYPDTWSTPHSYFSLTFCVQVQGKNKREKCLFMDQTS AKVTCHKDANVRVQARDRYSSFWSEWASVSCS
Purity	98%
Source	Yeast
MW	34.6 kDa
Entrez	281857
Expiration Date	6 months from date of receipt.