

Product Datasheet

Bovine IL-36RA Recombinant Protein (orb1820895)

Description	The Bovine IL-36 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-1F5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MVLSGALCFRMKDAALKVLYLHDNQLLAGGLQAGKVIKGEESVVPNRS�DAKL SPVILGVHGGSQLSCGTGQEP TLKLEPVNIMELYHSAEKSKKFTFYRRDTGLTS SFESAAYPGWFLCTVPEADQPLQITQLPKDTSWDNPIIDFYFQQCD
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	518514
Expiration Date	6 months from date of receipt.