

Product Datasheet

Swine CCL3L1 Recombinant Protein (orb1820899)

| | |
|-------------------------|---|
| Description | The Swine CCL3L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | CCL3L1 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | APLGADTPTACCFSYTSRQLPRKFVADYFETSSQCSPGVIFQTKKGREVCANPE DAWVQEYISDLELNA |
| Purity | 98% |
| Source | Yeast |
| MW | 7.8 kDa |
| Entrez | 494459 |
| Expiration Date | 6 months from date of receipt. |