

Product Datasheet

Swine CCL4 Recombinant Protein (orb1820910)

Description	The Swine CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTSCCFYTVRKLPRNFVTDYYETSSLCSQPAVVFQTKKGRQVCANP SDDWVQEYMDDLELN
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	396668
Expiration Date	6 months from date of receipt.