

Product Datasheet

Equine CCL3 Recombinant Protein (orb1820913)

Description	The Equine CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VPGADTPTACCFSYVSRQIPRKFINDDYETSSQCSKPAIIFQTKRSRQVCADPSE AWVQEYVTDLELSA
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	100057909
Expiration Date	6 months from date of receipt.