

Product Datasheet

Bovine IL-17A Recombinant Protein (orb1820932)

Description	The Bovine IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KAGVIIPQSPGCPPTEDKNFPQHVRVNLNIVNRSTNSRRPTDYHKRSTSPWTLHR NEDPERYPSVIWEAKCSHSGCINAEGKVDHMHNSVTIQQEILVLRRESQHCPS FRLEKMLVAVGCTCVTPIVRHLA
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	282863
Expiration Date	6 months from date of receipt.