

Product Datasheet

Bovine TNF alpha Recombinant Protein (orb1820933)

Description	The Bovine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRSSSQASSN KPVAHVADI NSPGQLRWWD SYANALMANG VKLEDNQLVV PADGLYLIYS QVLFRGQGCP STPLFLTHTI SRIAVSYQTK VNILSAIKSP CHRETPWAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPD YLDYAESGQV YFGIIAL
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	280943
Expiration Date	6 months from date of receipt.