

## Product Datasheet

### Equine IL-15 Recombinant Protein (orb1820942)

<b>Description</b>	The Equine IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	IL-15
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	NWQDVISDLK RIEDLIQSIH VDATLYTESD AHPNCKVTAM KCFLLELHVI SHESRNEDIK ETVENLIILA NSSLSSNGNV TESGCKECEE LEEKNIKEFL QSFVHIVQMF INPS (114) [Leu51Ser]
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	13.0 kDa
<b>Entrez</b>	<b>100034058</b>
<b>Expiration Date</b>	6 months from date of receipt.