

Product Datasheet

Bovine IL-2 Recombinant Protein (orb1820948)

Description	The Bovine IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTGNTMKEVKSLLLDLQLLLEKVKNPENLKLSRMHTFDYVPKVNATELK HLKCLLEELKLEEVNLNAPSKNLNPREIKDSMDNIKRIVLELQGSETRFTCEYDDA TVNAVEFLNKWITFCQSIYSTMT
Purity	98%
Source	Yeast
MW	15.5 kDa
Entrez	280822
Expiration Date	6 months from date of receipt.