

Product Datasheet

Chicken IL-16 Recombinant Protein (orb1820961)

Description	The Chicken IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SSTSSVASDA SQESTTEETI CTITLDKTAA GLGFSLEGGK GSIHGDKPII INRIFKGTSL EQSSPVQPGD ELLQVHTTAL QGLTRFEAWN IIKALPDGPI TAIKRKNPS SVTKKASETL
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	374270
Expiration Date	6 months from date of receipt.