

Product Datasheet

Bovine IL-6 Recombinant Protein (orb1820972)

Description	The Bovine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPLGEDFKND TTPGRLLLT PEKTEALIKR MVDKISAMRK EICEKNDECE SSKETLAENK LNLPMKEEKD GCFQSGFNQA ICLIRTTAGL LEYQIYLDYL QNEYEGNQEN VRDLRKNIRT LIQILKQKIA DLITTPATNT DLLEKMQSSN EWWKNAKIIIL ILRNLENFLQ FSLRAIRMK
Purity	98%
Source	Yeast
MW	20.7 kDa
Entrez	280826
Expiration Date	6 months from date of receipt.