

Product Datasheet

Swine IL-13 Recombinant Protein (orb1820990)

Description	The Swine IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LTCFGGLASP GPVPPHSTAL KELIEELVNI TQNQKTPLCN GSMVWSVNL TSMYCAALES LINISDCSAI QKTQRMLSAL CSHKPPSEQV PGKHIRDTKI EVAQFVKDLL KHLRMIFRHG
Purity	98%
Source	Yeast
MW	13.2 kDa
Entrez	396721
Expiration Date	6 months from date of receipt.