

Product Datasheet

Bovine IL-13 Recombinant Protein (orb1821002)

Description	The Bovine IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPVPSATALK ELIEELVNIT QNQKVPLCNG SMVWSLNLTS SMYCAALDSL ISISNCSVIQ RTKKMLNALC PHKPSAQVS SEYVRDTKIE VAQFLKDLLR HSRIVFRNER FN
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	281247
Expiration Date	6 months from date of receipt.