

Product Datasheet

Bovine IL-4 Recombinant Protein (orb1821004)

Description	The Bovine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GHKCDIT LAEIIKTLNI LTTRKNSCME LPVADVFAAP KNTTEKETFC RVGIELRRIY RSHTCLNKFL GGLDRNLNSL ASKTCSVNEA KTSTSTLKDL LERLKTIMKE KYSKC
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	280824
Expiration Date	6 months from date of receipt.