

Product Datasheet

LbuCas13a, Active, Glycerol Free (orb1821057)

| | |
|------------------------|--|
| Description | LbuCas13a, Active, Glycerol Free |
| Preservatives | Recombinant protein stored in 50mM Tris, pH 7.5, 600mM NaCl, 0.1mM TCEP. |
| Form/Appearance | liquid |
| Concentration | 0.1 µg/µl |
| Storage | Store at -70°C. |
| Tag | tag-free |
| Note | For research use only |
| Source | Leptotrichia buccalis |
| MW | 138.5 |
| Uniprot ID | C7NBY4 |
| Expiration Date | 6 months from date of receipt. |