

---

## Product Datasheet

**Lamin A+C Monoclonal Antibody [JOL2] -  
Fluorescein (FITC) / 50ul (orb1821243)**

**Description**

Lamin A/C Monoclonal Antibody:

<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Monkey
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, WB
<b>Immunogen</b>	Recombinant human lamin A
<b>Target</b>	This antibody reacts with both recombinant and native Lamin A and C in humans
<b>Preservatives</b>	Purified antibody conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% sodium azide
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	p02545
<b>Entrez</b>	<b>103223838, 4000</b>
<b>Dilution Range</b>	Optimal antibody dilution should be determined by titration, however as a guideline use at: WB 1:1000, IF 1:100, Flow cytometry 1:1000
<b>Expiration Date</b>	12 months from date of receipt.