

## Product Datasheet

# Extreme Thermolabile UNG (Uracil N-Glycosylase) (orb1821419)

<b>Description</b>	Extreme Thermolabile UNG (UDG) is ideal for use in one-step RT-PCR systems to remove contaminations ...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	1 unit/ $\mu$ l
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.