

## Product Datasheet

### Taq Polymerase - glycerol-free (orb1821432)

<b>Description</b>	Taq Polymerase-glycerol-free is recommended for use in freeze drying applications where glycerol mus...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	5 units/ $\mu$ l
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.