

## Product Datasheet

### Taq Polymerase (orb1821437)

<b>Description</b>	The polymerase is recommended for routine PCR applications (up to 4 kb fragment length) or high thro...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	5 units/ $\mu$ l
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.