

Product Datasheet

SYBR Green Fluorescent DNA Stain (orb1821443)

Description	SYBR Green Fluorescent DNA Stain is a superior DNA intercalator dye specially developed for DNA anal...
Form/Appearance	liquid, supplied in 20 mM Tris-HCl pH 8.5, 0.1 mM EDTA and 0.01 % Tween-20; Color: orange
Concentration	100 μ M
Storage	store at -20 °C. store dark
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 495 nm, λ_{em} 520 nm (bound to DNA).
Expiration Date	12 months from date of receipt.