

Product Datasheet

Estrogen receptor Antibody (orb1821539)

Description	Estrogen receptor Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The original antibody was generated by immunizing mice with estrogen receptor ER- α 36.
Target	Estrogen receptor
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The antibody specifically recognized ER- α 36 protein expressed by SUM159 cells, as confirmed by ELISA analysis (CN107987164).
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P03372
Expiration Date	12 months from date of receipt.