

## Product Datasheet

### Biotin Antibody (orb1821624)

<b>Description</b>	Biotin Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	The original antibody was generated by immunizing femal BALB/c mice with biotin-KLH conjugate.
<b>Target</b>	Biotin
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody was capable of binding biotin conjugated proteins in an ELISA. This antibody has a higher affinity for biotin conjugated protein than free biotin. This antibody is recommended for detection an diagnosis of biotin labeled proteins in a sample using techniques like flow cytometry, immunoblot assays, immunofluorescence studies.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.