

Product Datasheet

Biotin Antibody (orb1821626)

Description	Biotin Antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The original antibody was generated by immunizing femal BALB/c mice with biotin-KLH conjugate.
Target	Biotin
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody was capable of binding biotin conjugated proteins in an ELISA. This antibody has a higher affinity for biotin conjugated protein than free biotin. This antibody is recommended for detection an diagnosis of biotin labeled proteins in a sample using techniques like flow cytometry, immunoblot assays, immunofluorescence studies.
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.