

Product Datasheet

Cytomegalovirus (CMV) gH glycoprotein Antibody (orb1821667)

Description	Cytomegalovirus (CMV) gH glycoprotein Antibody
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	CMV antigen (Virus antigen batch V2547) prepared from MRC-5 cells infected with the Towne strain of CMV.
Target	Cytomegalovirus (CMV) gH glycoprotein
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q4WHS2
Expiration Date	12 months from date of receipt.