

## Product Datasheet

# Cytomegalovirus (CMV) gH glycoprotein Antibody (orb1821668)

<b>Description</b>	Cytomegalovirus (CMV) gH glycoprotein Antibody
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	CMV antigen (Virus antigen batch V2547) prepared from MRC-5 cells infected with the Towne strain of CMV.
<b>Target</b>	Cytomegalovirus (CMV) gH glycoprotein
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q4WHS2</a>
<b>Expiration Date</b>	12 months from date of receipt.