

## Product Datasheet

### STING Agonist-1 (orb1821676)

<b>Description</b>	Indirect activator of STING signaling
<b>Storage</b>	Storage Temp: -20°C
<b>Note</b>	For research use only
<b>Purity</b>	≥98%
<b>MW</b>	430.88
<b>CAS Number</b>	702662-50-8
<b>SMILES</b>	<chem>C1C(=O)N(C2=C(S1)C=CC(=C2)C(=O)NCC3=CC=CO3)CC4=C(C=CC=C4Cl)F</chem>
<b>Expiration Date</b>	12 months from date of receipt.