

Product Datasheet

HeV/NiV F/Fusion glycoprotein F0 Antibody [1F5] (orb1821720)

Description	HeV/NiV F/Fusion glycoprotein F0 Antibody [1F5]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Target	F
Preservatives	0.01M PBS, pH 7.4.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -80 °C for 12 months from the date of receipt.
Note	For research use only
Isotype	IgG1, kappa
Clonality	Recombinant
Purity	> 95% by SDS-PAGE.
Uniprot ID	Q9IH63 , O89342
NCBI	O89342 , Q9IH63
Expiration Date	12 months from date of receipt.