

## Product Datasheet

# HeV/NiV Glycoprotein G Antibody [HENV-26] (orb1821731)

<b>Description</b>	HeV/NiV Glycoprotein G Antibody [HENV-26]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, NeA
<b>Target</b>	G
<b>Preservatives</b>	0.01M PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -80 °C for 12 months from the date of receipt.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Uniprot ID</b>	<a href="#">O89343</a> , <a href="#">Q9IH62</a>
<b>NCBI</b>	<a href="#">O89343</a> , <a href="#">Q9IH62</a>
<b>Expiration Date</b>	12 months from date of receipt.