

Product Datasheet

Nipah virus/NiV Fusion glycoprotein/F Trimer Antibody [Fab92] (orb1821737)

Description	Nipah virus/NiV Fusion glycoprotein/F Trimer Antibody [Fab92]
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Target	F
Preservatives	0.01M PBS, pH 7.4.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -80 °C for 12 months from the date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Recombinant
Purity	> 95% by SDS-PAGE.
Uniprot ID	Q9IH63
NCBI	Q9IH63
Expiration Date	12 months from date of receipt.