

Product Datasheet

Nipah virus/NiV Prefusion Antibody [4H3] (orb1821744)

Description	Nipah virus/NiV Prefusion Antibody [4H3]
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Target	Nipah virus
Preservatives	0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 2 to 8 °C for 1-2 weeks. Store at -80 °C for 12 months from the date of receipt
Note	For research use only
Isotype	IgG1
Clonality	Recombinant
Purity	Protein A or G purified from cell culture supernatant - > 95% by SDS-PAGE
Uniprot ID	Q9IH63
NCBI	Q9IH63
Expiration Date	12 months from date of receipt.