

Product Datasheet

Nipah virus/HeV F/Fusion glycoprotein F0 Antibody (orb1821745)

Description	Nipah virus/HeV F/Fusion glycoprotein F0 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E. coli - derived recombinant Nipah virus/HeV F/Fusion glycoprotein F0 (Gly131-Thr546).
Preservatives	Unconjugated
Form/Appearance	Unconjugated
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -80 °C for 12 months from the date of receipt.
Note	For research use only
Application notes	ELISA:1:4000~1:8000, WB:1:5000~1:20000
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.