

Product Datasheet

BiTE Format Isotype Control for Bispecific Antibody (orb1821754)

Description	BiTE Format Isotype Control for Bispecific Antibody
Species/Host	Human
Conjugation	Unconjugated
Preservatives	Supplied in 0.01M PBS, pH 7.4.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store vial at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Application notes	Anti-Human CD3 & Anti-Lysozyme (HyHEL-10) Isotype Control
Isotype	VH-VL-VH-VL
Clonality	Monoclonal
Purity	> 95% by SDS-PAGE.
Source	XtenCHO
Expiration Date	12 months from date of receipt.