
Product Datasheet

Human EBI3 Protein, His Tag (orb1821761)

Description

Human EB13 Protein, His Tag is expressed

Lys 229 (Access...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

EB13

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.3 kDa. The protein migrates as 25-26 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

90%

Source

Human EB13 Protein, His Tag is expressed from E. coli cells. It contains AA Arg 21 - Lys 229 (Accession # Q14213).

MW

25.3 kDa

Uniprot ID
[Q14213](#)
Expiration Date

6 months from date of receipt.