

Product Datasheet

Anti-Influenza A [H1N1] Antibody (orb1821783)

Description	Anti-Influenza A [H1N1] Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A/California/14/2009 (H1N1) virions
Preservatives	PBS + 0.02% NaN3
Form/Appearance	Antigen Affinity Purified
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.