

## Product Datasheet

### Anti-PSDR1 (RDH11) Antibody (orb1821861)

<b>Description</b>	Anti-PSDR1 (RDH11) Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids from the N-terminal region of human PSDR1 protein
<b>Preservatives</b>	PBS + 0.02% NaN3
<b>Form/Appearance</b>	Antigen Affinity Purified
<b>Storage</b>	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q8TC12</a>
<b>Dilution Range</b>	WB: 1:500-1:1000
<b>Expiration Date</b>	12 months from date of receipt.