

Product Datasheet

Anti-Anti-His (C-terminal) Tag Antibody (orb1822153)

Description	Anti-Anti-His (C-terminal) Tag Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	Clone (M050) was generated from a recombinant protein containing a C-terminal His-tag.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1:5000, IHC: 1:100, ICC: 1:100, IP: 1:100
Expiration Date	12 months from date of receipt.