

Product Datasheet

Anti-Blasticidin S Deaminase Antibody (orb1822307)

Description	Anti-Blasticidin S Deaminase Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide (coupled to KLH) corresponding to amino acid residues in the N-terminal region of Blasticidin S Deaminase* from <i>Aspergillus Terreus</i> .
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13
Uniprot ID	POC2P0
Dilution Range	WB: 1:1000, ICC: 1:100
Expiration Date	12 months from date of receipt.